

FORA Advanced pro GD40h

SPECIFICATIONS



PRODUCT CODE	GD40h
CONNECTIVITY	Bluetooth V4.0
POWER SOURCE	2 x 1.5V AAA batteries
MEMORY CAPACITY	1,000 records
SYSTEM OPERATING CONDITION	10°C to 40°C
METER STORAGE CONDITION	-20°C to 60°C



PRODUCT CODE	ACS047
PIP CODE	409-3282
TEST STRIP	Blood Glucose Test Strips
ENZYME TYPE	GDH-FAD
MEASUREMENT RANGE	10 ~ 600 mg/dL (0.5 ~ 33.3mmol/L)
HEMATOCRIT RANGE	0 ~ 70%
CALIBRATION	Auto Coded
STRIP STORAGE CONDITION	2°C to 30°C
MAXIMUM SHELF LIFE (FROM DATE OF MANUFACTURE)	24 Months



PRODUCT CODE	ACS052
PIP CODE	409-3290
TEST STRIP	β-Ketone Test Strips
ENZYME TYPE	β-Hydroxybutyrate Dehydrogenase
MEASUREMENT RANGE	0.1 ~ 8.0 mmol/L
HEMATOCRIT RANGE	10 ~ 70%
CALIBRATION	Code Card
STRIP STORAGE CONDITION	2°C to 30°C
MAXIMUM SHELF LIFE (FROM DATE OF MANUFACTURE)	18 Months - Individually wrapped



FORA Advanced Pro GD40h

ACCURATE DUAL CONTROL

FORA Advanced Pro GD40h

ACCURATE DUAL CONTROL

FORA® Advanced pro series is a dual-parameter handheld device which measures both blood glucose and β -Ketone levels.

PROFESSIONAL DIABETES DEVICE WITH A RANGE OF FEATURES MAKES IT THE PERFECT PARTNER FOR ALL PEOPLE WITH TYPE 1 DIABETES:

<p>BLUETOOTH</p>  <p>CONNECTIVITY</p>	<p>WIDE HCT RANGE</p> <p>0-70% - Glucose 10-70% - Ketone</p>	 <p>FORA GD40h</p> <p>GLU AC β</p> <p>5.5 mmol/L</p> <p>6:59 AM</p> <p>M</p> <p>Exceeds Accuracy According to EN ISO 15197:2015</p>	<p>BACK LIT SCREEN</p> 	<p>STRIP FEED LIGHT</p> 	<p>DUAL PARAMETERS</p>  <p>AVERAGE RESULTS FOR 7, 14, 30, 60 & 90 DAYS</p> <p>■ The most cost effective β-Ketone test strips*</p>
<p>CONTROL SOLUTION</p>  <p>AUTO TAGGED</p>	<p>MAXIMUM</p>  <p>SHELF LIFE</p>		<p>SUITABLE FOR DRIVERS</p> 	<p>MEAL TAGGING</p> 	<p>FREE SELF-CARE APP</p>  <p>■ iFORA HM ■ BOLUS CALCULATOR</p>
<p>ISO</p> <p>EN ISO 15197:2015</p>	<p>MEMORY</p>  <p>1000 TEST</p>		<p>ONE-BUTTON DESIGN</p> 	<p>STRIP EJECTION</p> 	<p>SMALL SAMPLE VOLUME</p>  <p>■ 0.8μL for Blood Glucose Test ■ 1.0μL for β-Ketone Test</p>
					<p>QUICK RESULTS</p>  <p>■ 5 seconds for Blood Glucose Test ■ 10 seconds for β-Ketone Test</p>

* October 2019 Drug Tariff